



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FAMED ŻYWIEC LD-04 children's hospital crib (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Polish production,

Mattress with medical certificate in blue color,

Technical specifications:

Manual gradual tilt adjustment in the range: 0° / 13°,

Lowerable and adjustable side rails on both sides,

4 built-in bumpers,

Front / rear handles designed for free transport,

Mobile crib on 4 casters (2 braking casters),

Solid bed frame powder coated in white,

Includes:

Mattress,

Dimensions:

Overall: 82 x 150 x 138 cm,

Mattress: 68 x 124 x 9.5 cm,

Maximum load capacity: 40 KG,

Weight: 60 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

