



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET operator seat (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
German production,
New upholstery with medical certificate in the selected color,
The multifunctional MAQUET operator's chair is suitable for surgery, anesthesiology and outpatient care departments,
The stability of the chair is ensured by a stainless steel base with 5 twin wheels for hard floors,
Manual height adjustment in the range of 54 - 71 cm using a gas spring,
Manual backrest height adjustment in the range of 20 cm,
Manual backrest extension adjustment in the range of 7 cm,
Built-in foot stable central brake,
Dimensions:
Seat: 38 x 32 x 6 cm,
Backrest: 35 x 22 cm,
Weight: 32 KG,
It has a current inspection and is ready to work,
Technical Passport issued (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market,
Financing options (Poland only): Installments, Leasing, Loan,

