



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Micro-needle radiofrequency CELLINA RF DUAL CLINIC SYSTEM (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Korean production, (2021),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 300 VA,
The innovative CELLINA Radiofrequency creates a multifunctional platform, equipped with:
Non-invasive fractional matrix
64-pin RF MATRIX
LED

whose RF current energy produces a fractional effect, or pinpoint ablation of the skin, preserving the surrounding healthy, intact tissue. Heat-stimulated skin undergoes remodeling, followed by regeneration and thickening of connective tissue and rejuvenation,
Needle head with 25 automatic RF micro-needles -.

this advanced technology allows precise delivery of controlled RF energy directly to the skin through minimally invasive micro-needles. The gold-plated micro-needles automatically insert into the skin to a depth set by the operator and treat both superficial and deep layers of skin to achieve a 3D spatial effect. After the RF pulse, the needles automatically retract into the head,
Heads for vascular closure, treatment of active acne and enlarged pores,

Application:
Effective 3D lifting - thickening of all skin layers,
Smoothing wrinkles, stretch marks and scars,
Skin refreshing - brightening and even skin tone,
Treatment of hyperhidrosis,
Closing vessels, treating active acne and enlarged pores,

Features:
Tightens, tightens, thickens, and fills the skin with new collagen, creates a lifting effect,
Narrows pores,
Smooths wrinkles,
Brightens skin, evens skin tone,
Reduces redness and dilated vessels,
Shapes the jawline V - reduces the "double chin",
Remodels the cheeks,
Rejuvenates the eye area,
Rejuvenates the area around the mouth,
Smooths the texture of the epidermis,
Reduces acne,
Smooths the uneven surface of scars,
Reduces acne scars, including atrophic dimpled scars,
Reduces cellulite by tightening the skin layers,
Reduces stretch marks,

Specifications:
Display: TFT LCD,
Frequency: 1 MHz,
10 treatment levels,
Maximum power: 100 J/cm

2
,
Three technologies in one device,
Safe and bloodless treatment thanks to non-isolated needles,
Wide range of parameter adjustment allows treatments from very gentle to very intense,
Comfortable and lightweight treatment head allows for comfortable operation and visible treatment area,
Automatic RD needles with a thickness of 0.2 mm work evenly throughout the depth, enhancing the 3D effect of the treatment,
Precise injection depth setting (in 0.1 mm increments) with repetition setting from 0.5 to even 3 seconds,

Includes:
Power cord,
LED Micro Needle Head,
Matrix LED head,
PED head for vessel sealing,
Wired foot controller,
Dimensions: 43 x 72 x 100 cm,
Weight: 30 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

