



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ESAOTE MYLAB60 ultrasound machine (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Italian production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 1400 VA,
Stationary ultrasonic system
ESOATE MYLAB 60
is one of the more advanced ultrasound equipment. Particularly an advanced ultrasound system, diagnosticians love this machine because of its extremely customizable nature - allowing each doctor to configure the system according to their own preferences,
ESAOTE
designed the system with superior image enhancement technology, 3D/4D imaging, a wide range of potential imaging applications and compatibility with multiple probes, and then made sure the system was as user-friendly, configurable and cost-effective as possible,
Imaging technology available on the ultrasound system
MyLab
60 include VPan panoramic imaging, M-View composite imaging, Pure Brilliance (PBI) imaging and contrast-tuned imaging. Other features include X4D 4D volumetric imaging, trapezoidal imaging and real-time archiving,
Built-in DVD player,
Features:
Real-time archiving,
Integrated wireless connection,
Tp-View trapezoidal imaging,
Panoramic VPAN,
M-View (compound imaging),
PBI (pure brightness imaging),
CnTI (contrast tuned imaging),
X4D 4D volumetric imaging,
One-touch workflow,
50° rotatable control panel,
Height-adjustable control panel and LCD screens,
Compatible ultrasound probes/transducers:
MyLab CA631 curved array probe,
MyLab CA541 Curved Array AppleProbe,
MyLab CA431 probe with curved array,
MyLab CA430 probe with curved die,
MyLab CA1421 probe with curved die,
MyLab CA123 probe with curved die,
MyLab C5-2 R13 probe with curved die,
MyLab BC431 probe with curved die,
MyLab BL433 Linear AppleProbe array,
MyLab LA332 Linear AppleProbe array,
MyLab LA435 linear probe,
MyLab PA240 phased array probe,
MyLab PA121 probe with phased array,
MyLab BE1123 Endocavitary Probe,
MyLab EC1123 endocavitary probe
MyLab IOE323 Intraoperative probe in the shape of a field hockey stick,
MyLab TEE132 TEE probe,
MyLab Pencil CW 2 Pencil probe without imaging,
Menu in language: English,
Keypad on control panel in language: English,
Included:
ESAOTE LA523 linear head,
User manual in Polish,
Power cord,
Overall dimensions: 65 x 100 x 160 cm,
Weight: 90 KG,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

