



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



KLAROMED ceiling operating lamp plusLED 96/plusLED 56 (Double) (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production, (2014),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power of the lamp head plusLED 56: 26 W,
Power of plusLED 96 lamp head: 35 W,
Shadowless operating lamp
plusLED 96 / 56
Made with LED technology, lamp bowl in the shape of a circle. The body of the canopy monolithic, homogeneous, closed and compact, without segmentation gaps,
Canopy made of lightweight aluminum alloys for proper heat dissipation. Streamlined shape and special dome design for low aerodynamic drag. Canopy shape adapted to work with laminar air supply,
The light source is a white LED matrix, made up of only white LEDs with warm and cold color temperature. The LEDs are placed evenly inside the canopy in modules with the same properties, consisting of two LEDs and reflectors, i.e. mirrors that reflect the incident light from the LEDs, ensuring uniformity of illumination,
Technical specifications:
plusLED 96 PREMIUM:
Number of LEDs: 78,
Light fixture diameter: 70.8 cm,
Color temperature adjustable from 3500 K to 5500 K in 5 steps,
Light intensity at 1m distance: 160,000 lux,
Light intensity adjustable from 25÷100% ,
Operating field diameter: 160 - 360 mm,
Color reproduction index Ra: 99% ,
R9 color reproduction index: 99% ,
Depth of illumination (L1 + L2) at 20% : 175 cm,
plusLED 56 PREMIUM:
Number of LEDs: 54,
Light fixture diameter: 58.5 cm,
Color temperature adjustable from 3500 K to 5500 K in 5 steps,
Light intensity at 1m distance: 130,000 lux,
Light intensity adjustable from 25÷100% ,
Operating field diameter: 140 - 320 mm,
Color reproduction index Ra: 99% ,
R9 color reproduction index: 99% ,
Depth of illumination (L1 + L2) at 20% : 175 cm,
Technical parameters:
Grille width: 60 cm,
Mounting column height: 90 cm,
Length of freely adjustable left arm: 180 cm,
Length of freely adjustable right arm: 165 cm,
Included:
Mounting column,
Power transformer,
2 plusLED 56 and plusLED 96 lamp bowls,
2 mounting arms,
Grille,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,
The warranty period does not include bulbs and other light sources installed in the device,
Possibility of financing (Poland only): Installments, Leasing, Loan,

