



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Mobile surgical lamp COLD LIGHT YD01-5 (LED) (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
P.R.C. production, (2021),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 20 VA,
YD01-5 lamp
is designed for medium and small surgical operations. Easy to operate with free canopy height adjustment from 106 d0 174 cm,
Free intensity adjustment with built-in knob,
Adjustable with a knob to tilt the bulbs in the canopy,
Freely adjustable rotation of the canopy and its tilt,
Mobile on 4 casters (2 with brake),
Technical specifications:
Maximum light intensity: 80,000 lux,
Color temperature 5000±500K,
Depth of light intensity: 500 mm,
Total illumination 275 W/m²
,
Included:
Power cord,
2 pcs. - autoclavable canopy holders,
Dimensions:
Canopy diameter: ø49 cm,
Total: 48 x 100 x 195 cm,,
Total weight: 50 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,
The warranty period does not include bulbs and other light sources installed in the device,
Possibility of financing (Poland only): Installments, Leasing, Loan,

