



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Radiofrequency surgical diathermy ELLMAN PELLEVE S5 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Produced in the U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 120 W,
Pelleve S5 is a device that uses two types of radio waves for precise surgical procedures and skin lifting without a convalescence period,
Features:
Pelleve precisely delivers energy to the skin tissue using advanced radio frequency technology to induce collagen contraction without damaging the epidermis,
Compact ergonomic design allows for comfortable maneuvering around facial contours, including hard-to-reach areas such as the periorbital area,
Thanks to precise cutting, wounds heal much faster and are barely visible,
Can be used for various types of procedures,
Low risk of swelling after procedures,
Precise cutting in various tissues,
Digital displays of the set power,
Laminated, easy-to-use control panel,
Uses two types of radio waves for precise surgical procedures and skin lifting without a convalescence period:
Monopolar radio wave with a frequency of 4.0 Mhz,
Bipolar radio wave with a frequency of 1.7 Mhz,
It has the following functions:
Monopolar:
CUT - 120 W,
BLEND - 90 W,
COAG - 60 W,
FULGURATE - 45 W,
Bipolar:
CUT - 120 W,
BLEND - 120 W,
COAG - 120 W,
FULGURATE - 120 W,
Diathermy placed on a mobile ELLMAN trolley equipped with:
Convenient handle for transport,
Movable arm with clip,
Three shelves for accessories,
Four double wheels, including two front with brakes,
Included:
Handles:
Pelleve GlideSafe 20 mm. - 2 pcs. (1 new),
Pelleve GlideSafe 15 mm. - 2 pcs. (1 new),
Pelleve GlideSafe 10 mm. - 2 pcs. (1 new),
Pelleve GlideSafe 7.5 mm. - 3 pcs. (2 new),
New Ellman handle with switch - 3 buttons,
Ellman with switch - 3 buttons,
Ellman handle,
FLUKE 62 MINI laser thermometer,
ELLMAN flat electrode,
3-pin connecting cables - 4 pcs.,
Foot switch,
Power cord,
A set of Allen keys and accessories,
Overall dimensions: 50 x 35 x 130 cm,
Weight: 40 KG,
Cat 02,
It has a current inspection and is ready to work,
Technical Passport (Service Report) issued valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic market and international,
Financing options (Poland only): Installments, Leasing, Loan,

