



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **BODY-JET water liposuction (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
BODY-JET

is an innovative method to permanently remove unwanted fatty tissue throughout the body. It uses a modern method of water-jet assisted liposuction (WAL: water-jet assisted liposuction). The WAL technique used in the Body-Jet® device represents a breakthrough in previous liposuction procedures such as abdominal liposuction or knee liposuction. The focused jet of water separates fat tissue from fibrous tissue and muscle without significantly affecting the continuity of these structures.

The use of the WAL technique makes the Body-Jet procedure much less invasive than classic liposuction or even laser liposuction. The content of blood and other fluids in the aspirated fat is negligible, about 1-3% - in comparison, reaction fluids after classical liposuction are up to 40% . Using this innovative technique, we can also observe much less swelling, which usually accompanies this type of procedure. As a result, the patient can see the effect immediately after the liposuction procedure, and the small swelling makes it easier to precisely model the figure.

Selectable menu in the following language: English,

Keypad on the control panel in: English,

Included:

English instruction manual (PDF),

2 jars by Serres with a capacity of 3 liters each,

Foot attachment (2 pcs.),

Power cord,,

Dimensions: 53 x 55 x 179 cm,

Weight: 47 KG,

It has a current inspection and is ready to work,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,



