



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STEELCO DS 600D Disinfecting Washer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Italian production,
Power supply: 415 V,
Frequency: 50/60 Hz,
Power: 8250 W,
Circulating pump power: 750 W,
Dryer capacity up to 150 m
3
/h,
Dryer power: 4000 W,
Features:
Device for washing and thermodisinfection of surgical instruments, microsurgery instruments, anesthesia instruments and other medical devices,
Washer equipped with double-sided door designed for pass-through function,
The device entirely made of acid-resistant steel,
Washing chamber, washing arms and filters are made entirely of acid-resistant steel AISI 316L. The frame and housing are made of AISI 304 stainless steel,
Safe and easy-to-use control system with LCD display showing current cycle status,
Fully glazed door for visual inspection of the cleaning process,
The unit provides excellent cleaning and disinfection and drying results,
Injection of water and drying air through the same nozzles,
Forced air drying - with time and temperature control helps complete drying,
Ability to store up to 4 5-liter detergent canisters in the lower space of the unit,
Telescopic rails allow easy and safe loading/unloading,
Three-stage filtration system captures solid contaminants and prevents re-circulation,
Built-in water softener, increases the efficiency of the washing process,
Technical specifications:
Chamber capacity: 200 liters,
Garbage can capacity: 170 liters,
Unit with capacity of 6 to 8 DIN trays,
Thermal disinfection at 93°C and thermo-chemical disinfection,
Temperature is monitored by two independent sensors,
Tilt-open chamber door - when opened, it forms a work surface at a height of 850 mm,
Three basic washing programs: fast, standard and intensive,
Ability to modify programs (time, temperature, amount of detergent dispensed) to individual requirements,
Built-in ST2 printer,
1 detergent pump,
Noise level: 56 dB(A),
Included:
Power cord (415 V),
Batching cart having 3 shelves,
4 strainer baskets in size: 23.5 x 23.5 x 5.5 cm,
2 strainer basket with lid for small items in size: 10 x 10 x 4 cm,
2 pcs. - 8 position tool tray rack (4 cm spacing) in size: 18 x 39 x 11 cm,
8 instrument sterilization cassette (10 places): 18.5 x 14 x 1.5 cm,
8 instrument sterilization cassette (5 places): 18.5 x 8 x 1.5 cm,
Dimensions: 65 x 71 x 170 cm,
Weight: 225 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

