



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ASSISTINA 301 PLUS dental handpiece lubricator (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Austrian manufacture, (2019),  
Perfect maintenance in one move:  
The maintenance process starts automatically at the press of a button and ends after 35 seconds  
Automatic lubrication of the mechanical components of the tip that are in motion ensures optimal oil distribution, leaving an oil film,  
Contaminants are separated from metal components and removed. The life of the instrument is significantly increased and service costs are reduced,  
Product advantages:  
Cleaning and lubrication of internal mechanical components in motion during the cycle with service oil,  
Perfect cleaning of spray channels and air channels with cleaning fluid,  
Blowing with compressed air to remove excess oil,  
Includes:  
Instruction manual in Polish (PDF),  
Hose with cubic coupling,  
Adapter,  
Nozzle,  
Dimensions: 19 x 41 x 22 cm,  
Weight: 2.7 KG,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

