



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



METRUM CRYOFLEX CRYO-S CLASSIC cryosurgery apparatus (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Polish production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- The METRUM CRYOFLEX CRYO-S CLASSIC device is a model of a cryosurgery device. It is used for fast and effective tissue freezing using spray and contact cryoprobes. It is a pneumatically controlled device - it does not require electrical power from the mains. The freezing agent is nitrous oxide (N2O) or carbon dioxide (CO2) in steel pressure cylinders,
- Technical specifications:
 - o Device operating pressure: 3.5 - 5 Mpa,
 - o Maximum pressure: 5.5 Mpa,
 - o Cooling temperature: -89 (°C),
 - o Maximum device operation time from one tank: 240 min.,
 - o Built-in electronic treatment timer,
 - o Smooth manual pressure regulation,
 - o Used, among others, in: Surgery, Gynaecology, Dermatology, Phlebology,
 - o Mobile attached to a cylinder trolley (trolley on four wheels, two of which with brakes),
- Included:
 - o Power cord,
 - o New SD-10 Probe,
 - o New pressure hoses,
 - o Dedicated trolley,
 - o 10 l steel pressure cylinder,
 - o Dedicated case with accessories,
- Dimensions: 40 x 50 x 125 cm,
- Weight: 19 KG,
- Has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - o 6 months for the domestic market (Poland),
 - o 3 months for the international market,
 - o Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,
- Financing option (Poland only): Installments, Leasing, Loan,

