



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MAQUET RADIUS gynecological chair 1557.00AJ (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 230 VA,  
Maquet Radius  
is used for gynecological and urological examinations and procedures,  
Upholstery with medical approval in blue color (it is possible to replace the upholstery with a new one from our color pattern at an additional cost),  
Possibility of electromechanical adjustment with foot control:  
Tilt angle:  
Back segment in the range (- 0°) / 0° / 57°,  
Seat in the range of 25°,  
When adjusting the angle of the back segment, the seat angle is automatically adjusted,  
Storage compartment for 38 cm wide roll of backing paper,  
Cutout in the upholstery of the back segment for backing paper,  
Pull-out waste bowl with Ø 21 cm,  
Ability to electrically unfold the chair into the reclining position,  
Ability to freely adjust the foot rest vertically and 45° horizontally,  
Built-in patient step,  
Included:  
Wired footrest remote control,  
Waste bowl,  
Power cord,  
Dimensions:  
Back segment: 91 x 61 cm,  
Seat: 61 x 31 cm,  
Length in extended position: 175 cm,  
Overall: 155 x 76 x 138/168 cm,  
Maximum load capacity: 180 KG,  
Weight: 94 KG,  
Has up-to-date inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

