



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VEVOR YM850-40L Dental Compressor (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
P.R.C. production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 850 W,
Air compressor
VEVOR YM850
has a 40-liter tank to maximize efficiency. In addition, it is equipped with an oil-free pump for minimal maintenance,
The 1.1-hp/850-watt copper motor is very efficient and quiet, and has a heat dissipation hole to prevent overheating,
The compressor is equipped with a solenoid valve, a pressure gauge, a water separator and a safety valve to regulate the air pressure in the tank and keep it within a safe range at all times,
The oil-free motor operates with a noise level of 58 dB, which is not disturbing and provides a pleasant and calming working atmosphere,
The compressor has a handle for easy carrying and transportation, 4 non-slip feet provide stability during operation,
Technical specifications:
Housing material: steel,
Power: 1.1 HP / 850 W,
Tank capacity: 40 L. / 8.8 gallon,
Lubrication type; oil-free,
Speed: 1440 rpm,
Maximum pressure: 8 bar,
Noise level: 56 dB,
Overall dimensions: 43 x 43 x 71 cm,
Weight: 31 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

