



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GPE PSYCHOTRONICS STE-03 Stereometer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Device for diagnostics of spatial differentiation (stereoscopy), as well as visual acuity,
The full range of output positions of the measuring rods allows multiple measurement of near and far vision in an automatic manner,
The course of the examination including the current position of the rods is visualized on the monitor of the control device,
Features:
Variables measured:
Deviation of the measuring rods from the center value,
Number of underestimated and overestimated measurements,
Statistics for the left and right eye,
The program installed in the device performs statistical pre-processing of the results,
Form of presentation of results:
Average results,
Partial results,
Graphical presentation in the form of graphs,
The device cooperates with the computer program STiAD,
Menu of the device in Polish language,
Included:
Power supply with power cord,
2x RJ-45 cables,
Control computer,
Control module,
Overall dimensions: 30 x 76 x 118 cm,
Weight: 10 KG,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

