



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDICART IMOTECH MEDICAL endoscopy cart (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
U.K. production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Specifications:
Pull-out top keyboard drawer,
Power switch on the front of the cart,
Equipped with built-in strip of 6 devices,
Handle on the back for free transportation of the cart,
Mobile on 4 casters (2 braking),
Included:
Power cord,
Dimensions:
Top shelf: 48 x 28 cm,
2x middle shelf: 44 x 48 cm,
Lower shelf: 44 x 52 cm,
Overall dimensions: 64 x 65 x 120 cm,
Weight: 25 KG,
It has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

