



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE VOLUSON E6 ultrasound machine (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Austrian production, (2016),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 800 VA,
Voluson E6

is a high-end ultrasound scanner with RealTime 4D 3D imaging. It is a system with measurement and calculation software applicable to all fields of research, especially obstetrics, gynecology and radiology,

Features:

Ergonomic design, electronically controlled desktop, 19" LCD flat screen,

HD Live,

VCI with OmnieView mode - a solution that optimizes the resolution and contrast of the acquired spatial images in any chosen plane.

Omniview mode allows selection of any plane of any shape in the obtained "volume",

Measurement automation: Sono NT, Sono IT, Sono Biometry,

STIC,

TUI (tomographic imaging), SRI, B-flow,

Anatomical M-mode,

Works with external 4DView and ViewPoint software,

Technical specifications:

Stationary axial,

Number of active ports for connecting imaging heads: 4,

THI harmonic imaging,

3D imaging,

4D imaging,

Doppler presentation: (PW/CW/CD/PD/PDD/HPRF),

Panoramic imaging,

DICOM software,

Color monitor,

Screen diagonal: 19" inches,

Possibility of archiving on internal disk,

Possibility to connect to a computer,

Built-in recorder,

Selectable menus in DU, NED, EN, FI, FR, DE, IT, JP, NR, ES, SE, CN,

Included:

Power cord,

SONY UPP-110 printer,

GE IC5-9-D endovaginal probe (frequency 4 - 9 MHz),

CONVEX GE RAB6-D head (frequency 2 - 9 MHz),

Dimensions: 60 x 98 x 156 cm,

Weight: 130 KG

It has an up-to-date inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for a period of 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

