



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **GE MD15T cardiac monitor with SOLAR 8000M modular monitoring system (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Monitoring system  
GE SOLAR 8000M  
with monitor  
GE MD15T  
serves as the center of the monitoring system, user controls include a keypad and control knob for easy navigation,  
Modular patient monitoring system enhances clinical decision-making with a continuous flow of high-quality information. It provides robust clinical standards, including basic vital signs and advanced cardiac care,  
Transmission of patient information is via the  
TRAM-RAC2,  
Technical specifications:  
TRAM-RAC 2 module:  
ECG,  
TEMP/CO,  
2x BP,  
BP/BP,  
NBP,  
DISPLAY  
Defib. Sync,  
SPO2,  
Wired modules:  
ECG - Factory new 3-lead ECG cable,  
SOP  
2  
- Factory new sensor,  
NIBP - Cuff + air line,  
Includes:  
Power cord,  
ECG - Factory new 3 lead ECG cable,  
SOP  
2  
- Factory new sensor,  
NIBP - Cuff + air line,  
Connection cable set,  
Power supply,  
Dimensions:  
Overall: 35 x 40 x 40 cm,  
Total weight: 20 KG,  
It has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

