



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Video image processor and light source OLYMPUS OTV-SI (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Main and emergency lighting: 2 halogen bulbs 15 V / 150 W,
Suitable for flexible and rigid endoscopes,
It has the following functions:
Operation modes: manual and automatic,
Contour enhancement facilitating tissue observation,
White balance,
Gradual adjustment of light intensity,
It has the following inputs and outputs:
2 x Y/C,
VIDEO OUT,
SYSTEM,
H. COPY,
VTR REMOTE,
Keyboard connection socket,
Keyboard on the control panel in: English,
Included:
OLYMPUS Keyboard Model MAJ-1124,
Cable power supply,
User manual in English (PDF),
Dimensions: 43 x 30 x 17 cm,
Weight: 13 KG,
Cat 05,
It has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,
The warranty period does not cover bulbs and other light sources installed in the device,
Financing options (Poland only): Installments, Leasing, Loan,

