



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CANON TX-F non-contact tonometer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 80 VA,
Technical specifications:
Built-in 5" monochrome CRT monitor,
Measurement range: 0-30/0-60mmHg,
Setting method: Automatic , semi-automatic and manual,
Measurement time: 0.003 seconds,
Working distance of 11mm to the cornea of the eye,
Fixation point: green LED light,
Has a built-in thermal printer (paper width 58 mm),
Gentler air blast for better patient comfort,
Joystick replaced by functional roller and trackball,
Electric adjustment:
Chin up 65 mm,
Canopy height up/down with trackball within 30 mm,
Left/right with trackball within 90 mm,
Front/rear using roller within 40 mm,
Built-in internal memory for storing up to 10 measurements of each eye,
Menu in language: English,
Included:
User manual in English (PDF),
Power cord,
Dimensions: 28 x 53 x 47 cm,
Weight:
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

