



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



QB-E2 Test Tube Block (Ex-Display)

New (post-show),
Realistic pictures of the product,
QB-E2 block
For microcentrifuge tubes 2.0 ml,
Depth: 35 mm., / width: 11 mm. spaces for test tube,
Dimensions: 16 x 5.64 x 7.11 cm,
Weight: 1 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

