



### **OLYMPUS TC-C2 endoscope cart (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
U.K. manufacture,  
Olympus TC-C2  
is a compact endoscope cart for transporting endoscopes and delicate medical equipment,  
Upper top with dimensions: 45 x 39 cm,  
Recessed lower top with dimensions: 45 x 37.5 x 6 cm,  
Pull-out keyboard shelf with dimensions: 38 x 37 cm,  
2 shelves measuring 45.5 x 39 cm,  
Monitor holder with infinite manual angle adjustment in a wide range,  
Endoscope holder with manual height adjustment within 14 cm,  
3 large handles for moving,  
Distance between upper and lower top: 100 cm,  
The edges of the cart over the wheels protected by rubber,  
Mobile on casters, including two with brakes,  
Included:  
Endoscope holder,  
Monitor holder,  
Power cord,  
Dimensions: 70 x 60 x 150 cm,  
Weight: 40 KG,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

