



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **GLOBAL MAXI 4000 ENT chair (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
Produced in the U.S.A., (2015),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Upholstery with a medical certificate in blue (possibility of replacing the upholstery with a new one, from our color sample for an additional fee),  
The GLOBAL MAXI 4000 armchair is a multifunctional armchair ideally suited to ophthalmological and laryngology offices, and thanks to the built-in functions it can be used as a treatment armchair,  
Technical specification:  
Electric control of the armchair through a built-in double-sided side panel:  
Backrest tilt in the range (automatic synchronization of the leg rest with the backrest): 80° / 0° / -8°,  
Height adjustment of the armchair in the range: 56 - 98 cm,  
Possibility of setting in a lying position,  
Built-in programmable button with the option of assigning a selected position,  
Built-in button with a programmed sitting position,  
Electric locking of the chair rotation around its own axis,  
Manual functions:  
Double-sided retractable armrests,  
Possibility of rotating the chair around its own axis in the range of 360°,  
Free adjustment of the patient's headrest,  
Included:  
Power cord,  
Dimensions:  
Headrest: 23 x 21 cm,  
Backrest: 58 x 62 cm,  
Seats: 56 x 44 cm,  
Legrests: 44 x 34 cm,  
Patient's step: 34 x 21 cm,  
Total: 80 x 75 x 110/153 cm,  
Maximum load capacity: 180 KG,  
Weight: 150 KG,  
It has a current inspection and is ready for work,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,  
Financing options (Poland only): Installments, Leasing, Loan

