



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Fotel laryngologiczny SMR MAXI SELECT 28000-07 (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
Produced in the U.S.A.,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Upholstery with a medical certificate in blue (possibility of replacing the upholstery with a new one, from our color sample for an additional fee),  
The SMR MAXI SELECT armchair is a multifunctional armchair that is ideal for ophthalmological and ENT offices, and thanks to the built-in functions it can be used as a treatment chair,  
Technical specifications:  
Electric control of the armchair through a built-in double-sided side panel:  
Backrest tilt in the range (automatic synchronization of the leg rest with the backrest): 80° / 0° / -8°,  
Height adjustment of the armchair in the range: 53 – 96 cm,  
Possibility of setting the armchair in a lying position,  
Built-in programmable button with the option of assigning the selected position,  
Built-in button with programmed sitting position,  
Electric lock of the chair rotation around its own axis,  
Manual functions:  
Double-sided retractable armrests,  
Possibility of rotating the chair around its own axis in the range of 270°,  
Free adjustment of the patient's headrest,  
Possibility of retrofitting a treatment lamp on the back of the chair (retrofitting the lamp for an additional fee),  
Included:  
Power cord,  
Dimensions:  
Headrest: 25 x 23 cm,  
Backrest: 48 x 60 cm,  
Seats: 50 x 45 cm,  
Legrests: 42 x 32 cm,  
Patient's step: 28 x 21 cm,  
Total: 80 x 75 x 110/153 cm,  
Weight: 120 KG,  
It has a current inspection and is ready to work,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,  
The warranty period does not cover bulbs and other light sources installed in the device,  
Financing options (Poland only): Installments, Leasing, Loan,

