



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OASIS BATH WARMER infusion fluid warmer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
U.K. production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 200 W,
Infinitely adjustable heating temperature (max 60 °C),
Removable insert with five compartments (6.5 x 27.5 x 13 cm),
Includes:
Insert with five compartments,
Dimensions:
Compartments: 35.8 x 29.5 x 18 cm,
General: 41 x 36.5 x 29 cm,
Weight: 7.5 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan

