



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER VAPOR 2000 Inhalation Vaporizer (Reconditioned)

Product temporarily unavailable

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

German manufacture, (2012),

Drager Vapor 2000

Is a safe and precise vaporizer of inhalants,

Product features:

The volume of the vaporizer is 300 ml, so that the entire contents of a standard 250 ml anesthetic bottle can be poured into the vaporizer, even when there is anesthetic left inside from the previous refill,

The vaporizer does not require calibration during use,

Unique transport function - vaporizers can be carried in any position, as the transport safety device prevents spillage of the anesthetic, so it can be used immediately after transport,

All filling systems have been designed specifically for specific inhalation anesthetics, which means the risk of mistakenly using the wrong anesthetic is eliminated,

The particular vaporizer uses the Drager DW - 2000 attachment,

Fill type: Drager Fill,

Adjustable sevoflurane concentration in the range of 0 - 8% ,

Included:

User manual in English (PDF),

Dimensions: 11 x 25 x 20 cm,

Weight: 9 KG,

Cat 02,

Has an up-to-date inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for a period of 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

