



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LABCOLD VACCINE portable laboratory refrigerator RPDF0012F (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product
U.K. production,
Power supply: 230 V / 12 V (car charger),
Frequency: 50/60 Hz,
Power: 95 W,
Portable laboratory refrigerator designed for transporting vaccines and medicines,
Easy to read digital temperature display with adjustable temperature so you can be sure that the temperature displayed corresponds to the temperature in the container,
Technical specifications:
Built-in display,
Adjustable temperature range: +2 °C to +10 °C,
Capacity: 12 l,
Please note that this product is not suitable for transporting blood or blood products,
Included:
Power cord,
Car power cord,
Safety insert,
Dimensions: 20 x 33.5 x 16.5 cm,
Weight: 4.5 KG,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan

