



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDICOM colposcope VC-102 (video track) (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Polish production,

Power supply: 230 V,

Frequency: 50/60 Hz,

LED lighting,

Ergonomic head design

MEDICOM VC-102

together with the straight-arm and mobile base guarantee easy and convenient operation,

Integral glare reduction system completely eliminates reflections arising on wet surfaces,

Integrated monitor

NEOVO LCD

displays the examination "live",

Technical specifications:

MEDICOM VC-102 head:

Freely adjustable magnification: 4x to 50x,

Working distance range: 10 - 30 cm,

Field of view range: 6 - 75 mm,

Dual-zone ring illuminator composed of ultra-bright LEDs,

Image retention function,

Autofocus function or manual adjustment,

RED FREE function - electronic infrared or negative filter,

SCALE function - magnification marker or object size,

Glare reduction,

Possibility to connect ultrasound,

Green negative filter,

Working range of the arm: 60 cm,

NEOVO X-19AV monitor:

Screen diagonal: 19" inches,

LCD technology,

Resolution: 1280 x 1024 px,

TN matrix,

Brightness: 250 cd/m

2

,

Dynamic contrast ratio: 3000:1,

Response time: 3 ms,

Includes:

Power cord,

Overall dimensions: 88 x 70 x 134 cm,

Weight: 45 KG,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,

The warranty period does not include bulbs and other light sources installed in the device,

Possibility of financing (Poland only): Installments, Leasing, Loan

