



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EKOM DK 50 DM dental compressor (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Slovak production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Technical data:
Capacity: 5 l.,
Output pressure: 3.5 bar,
Operation mode: continuous S1,
Oil-free,
Built-in audible and visual alarms for cooling failure,
2 output connectors,
Volume:
≤ 50 dB [A],
Included:
Power cord,
Supply hose,
Dimensions: 40 x 51 x 72 cm,
Weight: 55 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan

