



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS WM-NP1 endoscope cart (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Made in the U.K.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Olympus WM-NP1 is a stable and durable trolley adapted to carry endoscopes and delicate medical equipment,
It has:
Upper tabletop with dimensions: 47 x 38 cm,
Lower tabletop with dimensions: 47.5 x 56 cm,
2 shelves with dimensions 40 x 54 cm with manual height adjustment (every 1.5 - 2 cm) using Allen screws,
Handles for easy moving of the trolley,
Endoscope holder with manual height adjustment,
Monitor holder with height adjustment in the range of 20 cm (6 points every 3.5 cm),
Electrical installation, 12 power sockets for the transported medical equipment,
Product features:
Power switch on the front of the trolley,
Detachable rear wall of the trolley makes it easy to manage cables,
Possibility to hide cables in panels on the sides of the trolley,
Shelf corners protected with protective covers, and edges above the wheels with tires,
Mobile on four wheels with Ø 12 cm, including two with brakes,
Included:
Built-in power cable,
Overall dimensions: 67 x 66 166/200 cm,
Weight: 86.5 KG,
Cat 02,
It has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic market and international,
Financing options (Poland only): Installments, Leasing, Loan

