



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



GE DASH 5000 Patient Monitor/Cardiac Monitor (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Manufactured in the U.S.A.,
Power supply: 230 V,
Backup power supply: via two built-in 11.1 V / 3.5 Ah batteries (not included, can be fitted with new ones for an additional fee),
Frequency: 50/60 Hz,
Power: 75 W,
Specification:
Colour LCD display with a diagonal of 12.1" and a resolution of: 800 x 600 px., provides an efficient screen image from a long distance,
Sockets: Defib Sync. / LAN / Aux,
Custom display options allow:
View and analyze data in graphical and numerical formats,
Real-time compilation of measurements and trends,
Free configuration of every screen element - from waveforms to parameter names,
Fully flexible alarm settings enable personalized monitoring,
Bright bedside alarm light helps to immediately locate the patient in critical situations,
Both basic and low-flow sidestream CO2 measurement options for neonatal to adult patients,
Lightweight, durable and ergonomic design, optimized for demanding transport,
Features:
TRIM KNOB operation makes it easy to enter information and navigate the screen,
Easy-to-use interface,
Equipped with:
NIBP - New dual-airway NIBP cuff,
ECG - 10-lead ECG cable,
SPO2 - Clip sensor,
TEMP - New sensor temperature (not included - can be retrofitted for an additional fee),
PI1 and PI2 - Can be retrofitted with a measuring line for an additional fee,
Possibility to select the menu in the following languages: EN, FR, IT, ES, DE, NL, SE, DK, HU, NR, PA, PL, RU, CZ, JP,
Included:
User manual in English (PDF),
Power cord,
NIBP - New NIBP cuff with double air tube,
ECG - 10-lead ECG cable,
SPO2 - Clip sensor,
Dimensions: 22.5 x 31.5 x 30 cm,
Weight: 10 KG,
Cat 02,
It has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty: 6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,
Financing options (Poland only): Installments, Leasing, Loan,

