



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **GE DASH 5000 Patient Monitor/Cardiac Monitor (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
Manufactured in the U.S.A.,  
Power supply: 230 V,  
Backup power supply: via two built-in 11.1 V / 3.5 Ah batteries (not included, can be fitted with new ones for an additional fee),  
Frequency: 50/60 Hz,  
Power: 75 W,  
Specification:  
Colour LCD display with a diagonal of 12.1" and a resolution of: 800 x 600 px., provides an efficient screen image from a long distance,  
Sockets: Defib Sync. / LAN / Aux,  
Custom display options allow:  
View and analyze data in graphical and numerical formats,  
Real-time compilation of measurements and trends,  
Free configuration of every screen element - from waveforms to parameter names,  
Fully flexible alarm settings enable personalized monitoring,  
Bright bedside alarm light helps to immediately locate the patient in critical situations,  
Both basic and low-flow sidestream CO2 measurement options for neonatal to adult patients,  
Lightweight, durable and ergonomic design, optimized for demanding transport,  
Features:  
TRIM KNOB operation makes it easy to enter information and navigate the screen,  
Easy-to-use interface,  
Equipped with:  
NIBP - New dual-airway NIBP cuff,  
ECG - 10-lead ECG cable,  
SPO2 - Clip sensor,  
TEMP - New sensor temperature (not included - can be retrofitted for an additional fee),  
PI1 and PI2 - Can be retrofitted with a measuring line for an additional fee,  
Possibility to select the menu in the following languages: EN, FR, IT, ES, DE, NL, SE, DK, HU, NR, PA, PL, RU, CZ, JP,  
Included:  
User manual in English (PDF),  
Power cord,  
NIBP - New NIBP cuff with double air tube,  
ECG - 10-lead ECG cable,  
SPO2 - Clip sensor,  
Dimensions: 22.5 x 31.5 x 30 cm,  
Weight: 10 KG,  
Cat 02,  
It has a current inspection and is ready to work,  
Issued Technical Passport (Service Report) valid for a period of 12 months,  
Warranty: 6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,  
Financing options (Poland only): Installments, Leasing, Loan,

