



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **AMPLIVOX 2160 audiometer (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

U.K. production,

Power supply: 230 V,

Frequency: 50/60 Hz,

AMPLIVOX 2160

is a portable audiometer that is ideal for schools, industry, and primary medical practices or wherever quick hearing screening is performed,

Specifications:

Frequency range: 125 - 8000 Hz,

Volume range: -10 dB to 120 dB,

Masking:

Pure tone: 11 band tapers,

Speaking: full white tone range from -10 dB to 100 dB,

Interrupter:

Manual / continuous / pulsed,

Control panel keypad in language: English,

Included:

User manual in English,

Power cord,

Bone headphone,

In-ear headphone,

In-ear headphone,

Wired microphone,

Patient button,

Dedicated carrying case,

2 x key,

Overall dimensions: 36 x 36 x 16 cm,

Weight: 5 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

