



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



YONKER YK-8000C patient monitor/cardiac monitor (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
P.R.C. production, (2020),
Power supply: 230 V,
Emergency power supply: via built-in 11.1 V / 2200 mAh battery,
Frequency: 50/60 Hz,
Power: 70 VA,
Advanced multi-parameter cardiac monitor
YONKER YK-8000C
, which can be used during surgery, post-surgery recovery, post-anesthesia and rescue to monitor the vital signs of adults, children and newborns. YK-8000C is a handy - ergonomic, comfortable and intuitive to use device,
Cardiomonitor designed for adults, children and newborns. Age group selection in the settings with automatic adjustment of alarm limits,
Technical specifications:
Screen: 12.1" inch LCL TFT,
Built-in measurement modules:
ECG,
SPO2,
TEMP1 | TEMP2,
NIBP,
The cardiac monitor works with a central monitoring station,
Adjustable alarm limits,
Memory of measured values (tabular and graphical trends),
Included:
Power cord,
ECG 5-lead cable,
SPO2 sensor,
NIBP cuff with air tubing,
Dimensions: 30.5 x 16.2 x 29 cm,
Weight: 3.3 KG,
Cat 01,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

