



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Patient monitor/cardio monitor YONKER YK-8000C (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

Produced by P.R.C., (2019),

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 70 VA,

Advanced multi-parameter cardio monitor YONKER YK-8000C, which can be used during surgery, post-surgery recovery, after anesthesia and rescue to monitor vital signs of adults, children and newborns. YK-8000C is a handy - ergonomic, comfortable and intuitive to use device,

Cardio monitor intended for adults, children and newborns. Age group selection in settings with automatic adjustment of alarm limits,

Technical specifications:

Screen: 12.1" inch LCL TFT,

Built-in measurement modules:

ECG,

SPO2,

TEMP1 | TEMP2,

NIBP,

The cardiomonitor cooperates with the central monitoring station,

Alarm limit adjustment,

Memory of measured values (tabular and graphic trends),

Included:

Power cord,

5-lead ECG cable,

SPO2 sensor,

NIBP cuff with air tube,

Dimensions: 30.5 x 16.2 x 29 cm,

Weight: 3.3 KG,

Cat 02,

It has a current inspection and is ready to work,

Technical Passport (Service Report) issued valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,

Possibility of financing (Poland only): Installments, Leasing, Loan,

