



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Ultrasonic cleaner 10L BS-UC10 (Ex-Display)

New (post-show),
Real photos of the product,
Manufactured by P.R.C.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Professional ultrasonic cleaner
BS-UC10

. Ultrasonic washer for disinfection of cosmetic accessories, especially microdermabrasion heads, cutters, etc., as well as dental and cosmetic tools,

The washers have a wide range of applications and will prove themselves in cosmetology and medicine,

In addition to very high cleaning efficiency and reaching the most inaccessible details, ultrasonic washing causes absolutely no mechanical damage to the washed objects, thanks to the high frequency of cleaning and the technology of heating the cleaning liquid, it allows you to thoroughly remove even the smallest dirt,

Durable construction of the washer, is made of stainless steel,

Technical specifications:

Capacity: 10 l.,

Ultrasonic power: 240W,

Frequency : 40KHz,

Operation time: 1-30 minutes,

Heating power: 300W,

Temperature control: 0-80 °C,

Included:

Power cord,

Mounting basket,

Dimensions:

Chambers: 30 x 27 x 15 cm,

Overall: 33 x 27 x 31 cm,

Weight: 3 KG,

Has up-to-date inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

