



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MEDIMA S300 infusion (syringe) pump (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Review photos,

Polish production, (2022),

Power supply: 230 V,

Emergency power supply: through 1 built-in 9.6 V / 2.2 Ah batteries,

Frequency: 50/60 Hz,

Power: 20 VA,

Medima S300

has an expanded drug library with limits for each infusion parameter. Soft limits (soft limits, min max) warn when the infusion parameter exceeds the recommended values, but do not limit the possibility of making changes. Hard limits (hard limits, min max) block the possibility of entering values outside their range,

Product features:

Large 3.2" color graphic display with high resolution and contrast,

Intuitive user interface,

Touch screen and numeric keypad,

Device menu in Polish language,

Extensive alarm system (allows you to adjust the volume and sound mode, night mode with the possibility of setting a muted sound and reduced display brightness),

Automatic syringe installation and recognition system,

Wide range of syringes (volumes of 2, 2.5, 3, 5, 6, 10, 12, 20, 30, 50 ml),

Very short time to start infusion,

Exceptional uniformity of infusion,

Advanced and intuitive drug library:

One for the entire hospital, with a total capacity of up to 5,000 dosing procedures,

Can be divided into up to 40 branches (CCAs) containing up to 500 drugs and have up to 40 defined categories,

Soft and hard limits for individual parameters,

Wide selection of infusion types and infusion units (ng, µg, mg, µU, mU, U, kU, µmol, mmol, mol, kg / min / hr / 24 hr),

Automatic and manual bolus (ability to program dose, time or rate of delivery),

Automatic elimination of occlusion bolus,

Integration with hospital databases (HIS/PDMS),

Cooperation with specialized software modules such as Medima Drug Editor, Medima Configurator, Medima Service & Maintenance,

Specifications:

Flow rate:

0.1 - 99.99 ml/h with increments of 0.01 ml/h,

100 - 999.9 ml/h with increments of 0.1 ml/h,

1000 - 2000 ml/h with increments of 1 ml/h,

Volume and infusion time:

0.1 - 20000 ml; 1 min - 200 hrs,

Bolus, initial bolus:

Rate: 20 - 2000 ml/h,

Volume: 0.1 - 50 ml,

Automatic and manual bolus,

Occlusion:

12 levels (75 - 900 mmHg),

Included:

Power cord,

Dimensions: 21 x 36.5 x 11.5 cm,

Weight: 2.5 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

