



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS SHOCKPULSE SE SPL-G lithotripsy system (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
OLYMPUS SHOCKPULSE-SE
is a system that offers dual-action lithotripsy with a patented technology that enables simultaneous generation of constant ultrasound wave energy and intermittent shock wave energy (ballistic/mechanical) with a single probe.
ShockPulse
offers unparalleled performance and significantly reduces treatment time,
Frequency: 19,500 Hz - 21,500 Hz,
Features:
Higher speed and performance: Effectively shreds and sucks up stones of various sizes, shapes and compositions,
Intuitive control for the surgeon: Ergonomic handle with integrated suction and hand activation,
Versatility: Full range of probes and sizes (for kidney, ureter and bladder),
One-handed design, no need to connect pneumatic or ultrasound transducers,
Manual transducer activation eliminates the need for a foot switch,
Ergonomic design allows for streamlined control during PCNL, URS and bladder procedures,
Transducer maintains a lower surface temperature throughout the procedure,
Control panel keypad in English,
Included:
Power cord,
OLYMPUS SPL-T Handpiece,
OLYMPUS SPL-NC nose cone,
OLYMPUS EGSPL-FS Foot Controller,
OLYMPUS torque wrench,
Dimensions: 36 x 39 x 11 cm,
Weight: 7 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

