



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **AMPLIVOX 2150 diagnostic audiometer (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Danish production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
AMPLIVOX 2150 audiometer

is a portable audiometer that is ideal for schools, industry and primary care practices or anywhere quick hearing screening is performed,  
Easy to use, it requires minimal training prior to use and offers calibrated pure tone and tone stimuli over a wide frequency range,

Specifications:

Frequency range: 125 to 8,000 Hz,

Volume range: From -10 to 100 (dB),

Masking: From 0 to 100 (dB),

Modes:

Left side of the handset (L),

Right side of the handset (R),

Both sides (L + R),

Bone handset,

Free volume control,

Keypad on control panel in language: English,

Included:

Power cord,

In-ear headphones,

Bone earphone,

In-ear headphone,

Patient button,

Dedicated carrying case for transporting the instrumentation,

Instrumentation dimensions: 31 x 16 x 10 cm,

Total weight: 5 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

