



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDTRONIC LIFEPAK 20E Defibrillator (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Manufactured in the U.S.A., (2020),
Power supply: 230 V,
Backup power supply: via 1 built-in 10.8 V / 6.6 Ah battery - (battery not included - can be retrofitted for an additional fee),
Frequency: 50/60 Hz,
Power: 120 VA,
Lifepak 20E is a two-phase, portable defibrillator, which is a complete response system in acute cardiac conditions,
Technical specifications:
AED mode,
Synchronization,
ECG parameter monitoring,
Easy-to-read 6.2" display,
Built-in printer (paper width: 50 mm.),
Possibility to select the menu in the following languages: ENG, PL, GER, FR, ESP, SWE, IT, DUT, FIN, DAN, NOR, PT, BZ, CS, HUN,
Keyboard on the control panel in: English,
Included:
User manual in Polish (PDF),
Power cord,
Hard paddles,
3-lead ECG cable,
Dimensions: 26 x 24 x 21 cm,
Weight: 6.5 KG,
Cat D142.
It has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,
Possibility of financing (Poland only): Installments, Leasing, Loan,



