



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**Infrasonic liposuction N.I.L. System EVA SP 1 (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Belgian production,  
EVA SP 1  
is the basic system of the Eva sp range  
,  
created for those users who have a system for INFILTRATION (infiltration pump) and ASPIRATION (powerful surgical suction),  
The device regulates the supply of compressed air  
to the handle  
LIPOMATIC 3  
,  
thereby enabling  
infrasonic liposuction with nutation  
(NIL® technique),  
It has an infrasonic drive,  
Advantages of the device:  
Non-fat structures are respected,  
The duration of the intervention is shorter, as is the postoperative subsidence,  
The surgeon's work is easier and more precise,  
Improved access to difficult and delicate areas,  
The body undergoes less trauma (hematoma, bruising, reduced swelling),  
The patient experiences less pain,  
The system is easy to transport,  
The system is compatible with NLF System® 1 and 2 lipofilling equipment,  
Standard equipment:  
LIPOMATIC 3 head  
(can be retrofitted with the new EVAMATIC 5 head at an additional cost),  
1x Foot Attachment,  
5x Pressure Hose,  
5x Suction hose,  
6x Cannula:  
1x Mercedes 4x300,  
1x 3G Accelerator 3x250,  
1x 3G Accelerator 3x150,  
1x Viterbo 4x250,  
1x Infiltration 2.5x200,  
1x Infiltration 2.5x300,  
1x Cannula adjustment wrench,  
1x Electric power cord,  
1x Cleaning key,  
1x Lipomatic ® oil (100ml),  
Infrasonic Drive,  
Control system,  
Cannula set, Accessories,  
Safety system,  
For an additional fee, the device can be retrofitted with:  
Compressor DK50 2V Silent S DRY Monsun - with silencing enclosure and dehumidifier,  
Compressor DK50 2V DRY Monsun - without housing with dehumidifier,  
EUROMI NLF2 fat transfer kit (Fat Collector / Fat Transfer)  
EUROMI NLF1 fat transfer kit (Fat Collector / Fat Transfer)  
MEDELA Fat Transfer Kit (Fat Collector / Fat Transfer)  
Consumables,  
Specifications :  
Infrasonic drive,  
Easy transportability,  
Power supply:  
4A for 200 / 240V power supply,  
8A for 100 / 115V~ power supply,  
Included:  
Power cord,  
Foot attachment,  
Device dimensions: 42 x 36 x 21 cm,  
Weight: 4 KG,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

