



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **PULPMATIC UNO pulp pot macerator (Ex-Display)**

New (post-show),  
Real photos of the product,  
U.K. production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 1 KW,  
Dish Macerator  
PULPMATIC UNO  
fits perfectly into small interiors such as dirty rooms and isolation rooms. The unit is opened and closed without the use of hands, making loading a load safe and hygienic for staff operating the unit and patients,  
Coating  
MicrobeSafe+  
prevents cross-contamination because it is an antimicrobial coating, and the fast and sharp knives, in a precise and lightning-fast manner, grind the pulp into particles that drain freely into the drainage system,  
Benefits:  
Economical - the most economical macerator in its class. Low power and water consumption are the strengths of this device,  
Safe and hygienic - reduces cross contamination, thanks to built-in manual and foot sensors that automatically open and close the macerator lid.  
Efficient - working with the device increases the efficiency and productivity of workers, eliminating the need to clean pools and ducks,  
The versatile device fits perfectly into a small space, such as an isolation room, but also complements large wards,  
Features:  
Audio messages,  
Automatic flushing,  
Free waste drainage,  
Touchless operation,  
Various starting options,  
Automatic drain flushing,  
Specifications:  
1 to 2 dishes per cycle,  
65 second cycle,  
Volume: 59 dB,  
Capacity: 15 l,  
Reduced water consumption,  
Includes:  
Instruction manual in Polish,  
Power cord,  
Disinfecting liquid,  
Dimensions: 40 x 51 x 97.5 cm,  
Weight: 64 KG,  
It has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

