



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ERBE ICC 300-H Surgical Diathermy (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

German production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 300 W,

Infinitely adjustable power,

Digital displays of the set power,

Foiled, easy-to-use control panel,

Has the following functions:

Monopolar cutting - max. 300 W:

AUTO CUT - cutting with automatic voltage control,

HIGH CUT - cutting tissues with low electrical conductivity,

Four types of coagulation - max. 120 W:

SOFT COAG - precise control of the coagulation process and prevention of carbonization through the coagulation mode,

FORCED COAG - when fast and effective coagulation is needed,

SPRAY COAG - suitable for argon gas coagulation,

AUTO BIPOLAR - bipolar coagulation,

Included:

User manual in English (PDF),

Power cord,

Monopolar electrode cable with ERGO PEN BOWA holder - 2 buttons, stem 4 mm (no electrodes included, can be retrofitted at additional cost),

Cable with BOWA clip for disposable neutral electrodes

(no electrodes included, can be retrofitted at additional cost),

Bipolar cable, BOWA tweezers

(can be retrofitted with bipolar tweezers at additional cost),

Foot attachment,

Dimensions: 41 x 37 x 15.5 cm,

Weight: 12 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,



