



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Slit lamp with tonometer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
P.R.C. production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Halogen bulb: 12 V / 30 W,
Technical data:
Microscope:
Eyepiece with magnification: 12.5x,
Diopter adjustment: +3 / -5 D,
Free adjustment of the eyepiece spacing,
Slit illumination:
Slit width:
Slit diameter: 0.2 / 1 / 3 / 4 / 6 / 10 mm,
Slit rotation: +/- 90°,
Vertical adjustment: 20 mm.,
Light source: HALOGEN 12 V / 50 W,
Maximum light intensity: 350,000 Lux,
Free adjustment of light intensity using a knob,
Built-in filters:
Blue filter,
Red-free filter (green),
Neutral,
Bright,
Built-in fixation point,
Slit lamp equipped with KEELER 2414-P-5000 TONOMETER,
Slit lamp placed on a mobile HAAG-STREIT table with manual height adjustment of the top in the range of 74 - 90 cm,
Included:
Power cord,
KEELER 2414-P-5000 TONOMETER,
Overall dimensions: 66 x 43 x 140/170 cm,
Weight: 40 KG,
It has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for the period 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,
The warranty period does not cover bulbs and other light sources installed in the device,
Financing options (Poland only): Installments, Leasing, Loan,

