



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **TOM-PAG ELEGANT medical couch (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Polish production, (2013),  
Rehabilitation table  
TOM-PAG ELEGANT  
is designed for therapeutic procedures in hospitals, sanatoriums and other health care units,  
Technical specifications:  
Structure made of steel profiles covered with powder paint,  
Bed and headrest made of particle board covered with polyurethane foam upholstered with washable material,  
Two-section sleep surface,  
Infinite manual adjustment of the angle of the headrest,  
Adjustable in range:  
Headrest angle adjustment: - 60° / 0 / 45 °,  
Paper hanger,  
Upholstery with medical approval in blue color (it is possible to replace the upholstery with a new one from our color pattern for an additional fee),  
Dimensions:  
Headrest: 58 x 44 cm,  
Body segment: 58 x 144 cm,  
Overall length: 190 cm,  
Height: 52 cm,  
Maximum load capacity: 170 KG,  
Weight: 15 KG,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

