



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE DATEX-OHMEDA S5 AVANCE anesthesia machine (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Manufactured in the U.S.A.,
Power supply: 230 V,
Backup power supply: via 2 new built-in 12 V / 5 Ah batteries,
Frequency: 50/60 Hz,
Software: 3.30,
Datex-Ohmeda S/5 Avance is a device for anesthetizing neonates, children and adults with full monitoring of ventilation and concentration of medical gases and anesthetics,
It has two ports for vaporizers with the Selectatec (S-2000) mounting system - no vaporizer included (can be retrofitted for an additional fee),
Powered by gases from a central installation: O2, N2O, Air,
Possibility of supplying gases from additional cylinders: O2, N2O,
Electronically controlled gas mixer with pneumatic control of reserve ventilation,
Built-in color ventilation screen with a diagonal of: 12",
Built-in LED screen backlight with adjustable light intensity,
Easy access to the menu,
Intuitive user interface,
Full patient monitoring capabilities,
Advanced ABS system:
Minimum number of components and cable connections significantly reduces the risk of leaks or disconnections,
All components except the oxygen sensor and flow sensors are autoclavable,
ABS breathing system is latex-free,
ABS system belonging to the common platform of GE Datex-Ohmeda devices: Aisys, Avance, Aespire,
Ventilation options:
Volume control,
Pressure control,
PSVPro® (pressure support with apnea backup mode),
Synchronized intermittent mandatory ventilation (SIMV) (volume and pressure),
Electronic PEEP,
Tidal volume compensation,
Worktop dimensions: 45 x 32 cm,
Two drawers dimensions: 32.5 x 29 x 18 cm,
Possibility to select the menu in the following languages: PL, SWE, RUS, TUR, CS, GER, HUN, DUN, NOR, FIN, CHIN, DU, PT, JP, IT, ESP, GER, FR, ENG,
Device equipped with:
Patient monitor / cardiomonitor GE DATEX-OHMEDA D-LCC15-03:
Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Manufactured in the U.S.A.,
Power supply 230 V, and through a new built-in 6 V / 1.2 Ah battery,
Frequency: 50/60 Hz,
Designed for multi-parameter patient monitoring,
Custom display options allow viewing and analyzing data in graphical and numerical formats, compiling measurements and trends in real time,
Equipped with:
Display: 15" inches,
Type: Active color matrix TFT LCD,
Resolution: SVGA 800 x 600 px,
S/5 ComBar:
Anesthesia-specific command bar for intuitive access to monitoring functions,
Removable command bar allows the user to control the monitor nearby,
Graphical trends:
All parameters, selectable time from 30 minutes to 72 hours,
Numerical trends:
Continuously measured parameters sampled every 5 minutes, NIBP and PCWP values updated after each measurement. In addition, the marker triggers a numerical trend update,
Trend cursor:
Both in graphical and numerical trends,
Acoustic and visual alarms,
Possibility to select the menu in the following languages: ENG, GER,
M-CAIO anesthetic gas measurement module with D-FEND water trap,
DATEX ENGSTORM module equipped with:
M-NESTPR module equipped with:
ECG - 3-lead cable,
SPO2 - Clip sensor,
NIBP - Cuff (6 pcs.) + air tube,
T1 and T2 - 1 sensor in the set (possibility of retrofitting a second new one for an additional fee),
Included:
User manual in English (PDF),
3-lead ECG cable,
SPO2 clip sensor,
Patient temperature sensor,
NIBP cuff with air tube: 3 pcs. (adult cuff),
3 pcs. (children's cuff),
Gas supply hoses (N2O, AIR, O2),
Gas supply hoses with tips to be chosen by the customer (AGA or DIN),
Patient's respiratory system,
Sampling line for gas measurement,
Power cord,
Dimensions:
Total: 76 x 73.5 x 175.5 cm,
Drawers: 32.5 x 29 x 18 cm,
Total weight: 125 KG,
Cat 03,
On the domestic market (Poland), in the case of installation of the device within 30 days from the date of delivery of the equipment, the installation and training are free of charge, after exceeding this period these services will be paid,
It has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for the period 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market,
Financing options (Poland only): Installments, Leasing, Loan,

