



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rehabilitation - hospital bed FAMED ŻYWIEC LR-12 (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Polish production,

Power supply: 230 V,

Rehabilitation bed LR-12

is designed for adult hospital patients during diagnosis, monitoring and treatment of diseases,

The autoregressive function in the backrest segment reduces the pressure on the patient's body in the lumbar and thigh region,

Frequency: 50/60 Hz,

Power: 300 VA,

Advantages:

The bed provides the patient with a faster recovery, gives the patient a feeling of safety and comfort in treatment, facilitates rehabilitation,

Eliminates factors that hinder the work of medical personnel (e.g. moving heavy patients, overloading the spine),

The autoregressive function in the backrest segment reduces the pressure on the patient's body in the lumbar and thigh region,

It is coated with powder varnish resistant to UV radiation, mechanical damage and disinfectants,

It is easy to keep clean,

Has the ability to be placed in the cardiac chair position,

Allows you to control the electric functions of the bed even while rolling,

Provides full safety with a bed load of up to 170 kg,

Allows x-ray imaging,

It is equipped with a central wheel lock to guarantee easy and stable bed immobilization, as well as a wheel for straight travel to allow one person to maneuver the bed,

Technical specifications:

Electric adjustment by wired handheld remote control:

Height adjustment: 43 - 80.5 cm,

Backrest segment elevation angle: 70°,

Thigh segment elevation angle: 37°,

Lifting angle of shank segment: 20°,

Electric positioning of the bed in the cardiac chair position,

Manual CPR function of the backrest and Anti-Trendelenburg position of the bed,

Medical approved mattress in white,

Includes:

Power cord,

Drip stand,

Wired remote control,

White mattress with medical certificate

(can be replaced with a new one at an additional cost),

Dimensions: 96.5 x 216 x 43/80.5 cm,

Maximum load capacity: 185 KG,

Weight: 150 KG,

Cat 03,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

