



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PRAISTECH E0S E01 treatment table (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Upholstery with medical certificate in cream color (it is possible to replace the upholstery with a new one from our color pattern for an additional fee),
Height adjustment of the table is carried out with an electric actuator controlled by a handheld remote control,
Fully bearing height adjustment mechanism is used,
Individual parts of the table top are adjusted by gas springs,
Features and advantages of the table:
Gas spring from German company BANSBACH,
Electric actuator of German company DEWERT,
Upholstery of very good quality,
5 cm thick foam,
In the moving places of the table used steel ball bearings, which affects the unreliability of the table,
Frame elements made of the best Polish steel,
Technical parameters:
Electric height adjustment by wired handheld remote control in the range: 48 to 100 cm,
Manual adjustment of the headrest in the range: -70° / 0° / 40°,
Built-in side Euro rails (40 cm),
Included:
Power cord,
Wired remote control,
70 cm wide paper hanger,
Dimensions:
Headrest: 80 x 44 cm,
Bed: 82 x 152 cm,
Overall: 82 x 200 x 48/100 cm,
Maximum load capacity: 160 KG,
Weight: 88 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

