



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **KARL STORZ 26003ACA HOPKINS rigid endoscope (laparoscope) (Reconditioned)**

New, (exhibition),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

German production,

Technical specifications:

Diameter: 10 mm,

Viewing angle: 0°,

Autoclavable,

Built-in fiber optic light transmission,

Includes:

KARL SOTRZ autoclavable cassette - ENDOSKOPE,

Dimensions:

Working length: 31.5 cm,

Overall length: 39.5 cm,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for a period of 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

