



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **CONVEX SAMSUNG Medison C2-8 head (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Korean manufacture,

C2-8 probe

Has a frequency range of 2.00 to 8.00 MHz,

Application:

Abdominal,

Obstetric Gynecology,

Tissue structures,

Specifications:

Frequency: 2 - 8 MHz,

Transducer type: convex,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

