



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



SAMSUNG Medison LN5-12 linear probe (Ex-Display)

New (exhibition),
Real product photos,
Korean production,

The LN5-12 probe has a frequency range of 5.00 to 12.00 MHz,
Application:

Musculoskeletal system,
Vascular system,
Small body parts,

Specification:

Frequency: 5 - 12 MHz.,

Probe type: linear,

It has a current inspection and is ready to work,

Technical Passport (Service Report) issued valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market,

Financing options (Poland only): Installments, Leasing, Loan,

